



**HAL**  
open science

## The Diverse Family of Arp2/3 Complexes.

Javier Pizarro-Cerdá, Dror Shlomo Chorev, Benjamin Geiger, Pascale Cossart

► **To cite this version:**

Javier Pizarro-Cerdá, Dror Shlomo Chorev, Benjamin Geiger, Pascale Cossart. The Diverse Family of Arp2/3 Complexes.. Trends in Cell Biology, 2017, 27 (2), pp.93-100. 10.1016/j.tcb.2016.08.001 . pasteur-01457856

**HAL Id: pasteur-01457856**

**<https://hal-pasteur.archives-ouvertes.fr/pasteur-01457856>**

Submitted on 6 Feb 2017

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial - ShareAlike| 4.0 International License

## The diverse family of Arp2/3 complexes

Javier Pizarro-Cerdá<sup>1,2,3</sup>, Dror Shlomo Chorev<sup>4\*</sup>, Benjamin Geiger<sup>4</sup>  
and Pascale Cossart<sup>1,2,3</sup>

Institut Pasteur, Unité des Interactions Bactéries Cellules, Paris F-75015, France<sup>1</sup>;

INSERM, U604, Paris F-75015, France<sup>2</sup>;

INRA, USC2020, Paris F-75015, France<sup>3</sup>;

The Weizmann Institute of Science, Department of Molecular Cell Biology, Rehovot I-7610001, Israel<sup>4</sup>

\*Present Address: Department of Chemistry, University of Oxford, Physical and Theoretical Chemistry Laboratory, South Parks Road, Oxford OX1 3QZ, UK

### Correspondance:

Javier Pizarro-Cerda (javier.pizarro-cerda@pasteur.fr) and Pascale Cossart (pascale.cossart@pasteur.fr): Unité des Interactions Bactéries Cellules, Institut Pasteur, 25 rue du Docteur Roux, 75724 Paris Cedex 15 France. tel: +33(0)1 4061 3779, fax: +33(0)1 4568 8706

### Abstract

The Arp2/3 complex has so far been considered to be a single seven-subunit protein complex required for actin nucleation and actin filament polymerization in diverse critical cellular functions including phagocytosis, vesicular trafficking, lamellipodia extension and cytokinesis. The Arp2/3 complex is also exploited by bacterial pathogens and viruses during cellular infectious processes. Three recent studies suggest that some subunits of the complex are dispensable in specific cellular contexts, pointing to the existence of alternative Arp2/3 complexes containing other components such as vinculin or  $\alpha$ -actinin, as well as different isoforms or phosphorylation variants of canonical Arp2/3 subunits. This diversity should be considered when assigning specific Arp2/3 assemblies to different actin-dependent cellular processes.

## Introduction

The actin cytoskeleton is one of the main components of eukaryotic cells, not only providing the molecular basis for cellular morphogenesis and migration, but also participating dynamically in mechanical resistance to deformation, uptake of extracellular material, intracellular vesicular transport, cytokinesis and cell adhesion. The actin cytoskeleton also participates in the organization of complex cellular structures such as filopodia, lamellipodia and podosomes [1,2].

Polymerization of actin monomers into actin filaments requires the activity of cellular actin nucleators. The Arp2/3 complex, the first nucleator identified in eukaryotic cells, plays a central role in many cellular processes and is highly conserved from trypanosomes to the fission yeast and humans [3-5]. Other nucleators include formins, Spire, Cordon-bleu (COBL) and Leiomodins [6].

The Arp2/3 complex is composed of seven subunits [7] and it has been traditionally considered as a single entity, associated with the vast majority of cellular processes in which its function is required and has been studied. Three recent studies [8-10] in the same mammalian cell system reveal that diverse Arp2/3 complexes may regulate different cellular and pathogen-associated functions, raising the interesting possibility that Arp2/3 complex compositions may have been overlooked, paving the way for the identification of novel complexes associated to different actin polymerization-mediated processes.

## Discovery and functions of the classical Arp2/3 complex

The Arp2/3 complex was first isolated from *Acanthamoeba castellanii* during a search for ligands of the actin-binding protein profilin [11]. It contained seven proteins: the actin-related proteins Arp2 (44-kD) and Arp3 (47-kD) considered as 'unconventional actins', together with a 40-kD protein similar to a WD40  $\beta$ -propeller protein from *Dictyostelium discoideum*, and four additional proteins of 35-, 19-, 18- and 13-kD [11]. Subsequently, the Arp2/3 complex was also identified and associated with actin-rich structures in the fission yeast *Schizosaccharomyces pombe* [12] and in the budding yeast *Saccharomyces cerevisiae* [13]. In human cells, the Arp2/3 complex consists of Arp2 and Arp3, together with the Arp complex subunits ARPC1, ARPC2, ARPC3, ARPC4 and

ARPC5 [14]. While consensus exists concerning the Arp2 and Arp3 nomenclature, different names have been used in the literature concerning the other Arp2/3 complex subunits: a nomenclature proposal across species is presented in **Table 1**.

The function of the Arp2/3 complex was shown for the first time to be critical in triggering actin polymerization when it was isolated from a subcellular fraction of human platelets that sustained actin assembly by the bacterial pathogen *Listeria monocytogenes* [15]. The *L. monocytogenes* surface protein ActA activates the Arp2/3 complex to initiate actin polymerization, and was the first actin nucleation promoting factor (NPF) to be identified [16,17]. Several mammalian NPFs were subsequently identified, including WASP [18], N-WASP [19], Scar/WAVE [20], and cortactin [21] (see **Text Box 1**). The purified *A. castellanii* Arp2/3 complex was shown to nucleate the formation of actin filaments at 70° from other filaments [22]. A series of elegant microscopy and biochemistry investigations then definitively established the key role of Arp2/3 in actin polymerization and the formation of branched structures [23,24] (see **Text Box 2**).

As mentioned above, the function of the Arp2/3 complex is subverted by bacterial pathogens at different stages of their infectious processes [25]. The Gram-positive pathogen *L. monocytogenes* uses Arp2/3 not only to mediate intra- and inter cellular movements but also to trigger cellular invasion [26-28]. The Gram-negative pathogen *Shigella flexneri* also requires Arp2/3 function for actin-based motility [29] and for bacterial internalization within host cells [30]. Interestingly *S. flexneri* does not express an ActA-like protein but instead recruits on its surface, via the protein IcsA/VirG, the NPF N-WASP which in turn activates Arp2/3 to mediate actin-based motility [29,31]. The Gram-negative bacteria *Rickettsia parkerii* and *R. conorii* activate Arp2/3 during early stages of bacterial intracellular motility via a protein called RickA [32-34]. Moreover, *R. parkerii* requires Arp2/3 activity to invade diverse host cells [35]. Vaccinia virus is able to move at the surface of cells on actin-based structures [36], which requires the function of the Arp2/3 complex [37]. Other bacteria including Mycobacteria [38] and *Burkholderia thailandensis* also move via an actin-based motility requiring Arp2/3 functions [39-42].

In *S. pombe* and *S. cerevisiae*, deletions of genes encoding each of the subunits of the Arp2/3 complex cause severe growth defects or lethality [43,44], suggesting a major role for all subunits *in vivo*. In particular, Arp2/3 had been shown to be important for the formation and function of cortical actin patches where clathrin-mediated endocytosis takes

place [13]. Mammalian Arp2/3 complex was localized to regions of lamellipodial protrusion [14,45] and together with cofilin and other actin-binding proteins, was shown to control the organization and tread-milling of actin filaments in lamellipodia [23]. The Arp2/3 complex has been associated to other cellular functions requiring actin polymerization including phagocytosis [46], trafficking within and from the Golgi apparatus [47] as well as formation of focal adhesions [48]. The critical role of Arp2/3 in humans is highlighted by the Wiskott-Aldrich syndrome (WAS), a recessive X-linked genetic disorder characterized by mutations in the WAS protein (WASP), which is characterized by defects in the actin-rich immunological synapse between T cells and antigen presenting cells, leading to severe defects in immunological responses [49,50].

Initial detailed analysis of the contribution of each Arp2/3 complex subunit to actin polymerization, using *L. monocytogenes* ActA as a NPF in a baculovirus expression system in insect cells, indicated that only Arp2 and Arp3 are directly involved in actin polymerization, the role of the other subunits being less clear [51]. More recent structural evidence [52] confirms an initial prediction that ARPC2 and ARPC4 provide the main surface for interaction of the complex with the mother actin filament [51]. ARPC3 is proposed to form a bridge between Arp3 and the mother actin filament [52] but complexes lacking ARPC3 display minor functional defects [44,51]. While ARPC1 is supposed to make only minor contacts with the mother actin filament [52], complexes lacking this subunit are far less effective in actin nucleation, suggesting additional roles for ARPC1 including binding of NPFs [53]. ARPC5 was proposed to tether Arp2 to the rest of the complex [52].

Several reports also suggest a functional role played by phosphorylation of different subunits of the Arp2/3 complex. ARPC1 phosphorylation by p21-activated kinase (Pak1) was reported to be crucial for mammalian cell motility [54]. It has been suggested that Arp2 phosphorylation is required and critical for Arp2/3 complex binding to the pointed end of actin filaments and actin nucleation in cultured *Drosophila* cells [55], but mutation of the phosphorylated residues had only subtle effects on motility in *Dictyostelium* [56]. As shown recently, phosphorylation of Arp3 by the *Legionella pneumophila* kinase LegK2 inhibits actin polymerization at the surface of bacterial-containing phagosomes [57].

Several subunits of the Arp2/3 complex (i.e. Arp3, ARPC1 and ARPC5) display more than one isoform [14], but the functional significance of these variants had not been

investigated in detail. While the major subunit Arp3 is detected in all tissues, a gene encoding the isoform ARP3 $\beta$  was detected predominantly in brain neuronal cells and was proposed to play a role in the development and/or maintenance of nerve cells [58]. Two variants of ARPC1 presenting 70% homology had been known for long [12,45] and a mutation in the gene *ARPC1A* was shown to impact cell migration and invasion in pancreatic cancer [59]. ARPC5 was also found to display a second isoform, named ARPC5B, which exhibited a regular expression in many tissues but with the highest levels in the brain, while the original ARPC5A was found highly enriched in the spleen and thymus [60].

### **Diversity of Arp2/3 complexes**

**Focal adhesions.** Association of the Arp2/3 complex to focal adhesions in human skin cells had previously been shown to require interactions with vinculin [48]. A recent native mass spectrometry analysis of proteins extracted from the dense plaques (focal adhesion homologous structures) of chicken smooth muscle revealed surprisingly that Arp2/3 complexes present in these structures, as inferred from mass spectrometry results, are actually 'hybrid complexes', consisting of a core composed of Arp2, Arp3 and ARPC2, together with  $\alpha$ -actinin and vinculin, or Arp2, Arp3, ARPC2, ARPC3 and vinculin [8]. This study therefore supported, for the first time, the notion that alternative Arp2/3 complexes that do not consist of the seven classical subunits are involved in specific cellular processes. Notably, these alternative complexes contain vinculin that can mediate the recruitment of the complex to focal adhesions and compete with ARPC1B in HeLa cells; knock-down of ARPC1B has therefore a positive effect on focal adhesion and stress fiber formation, as the equilibrium is shifted towards Arp2/3-vinculin hybrid complexes formation [8].

***Listeria monocytogenes* infection.** Specific roles for ARPC1A and ARPC1B were recently identified in human genome-wide RNA interference (RNAi) screens investigating HeLa cell infection by *L. monocytogenes* [9]. Knock-down of ARPC1B but not of ARPC1A significantly diminished bacterial entry, highlighting a critical contribution for ARPC1B function in this context. Moreover, it was observed that ARPC4 and ARPC5 subunits do not contribute to *L. monocytogenes* cellular invasion. The contribution of the different Arp2/3 subunits to bacterial actin-tail formation was also studied, identifying a major contribution of Arp2, Arp3, ARPC1A, ARPC2, ARPC3 and ARPC4 in *L. monocytogenes*

actin-based motility, but no role for ARPC1B nor ARPC5 was found [9]. These results therefore show not only that ARPC5 is dispensable for both bacterial entry and actin-tail formation but also that ARPC1 isoforms contribute to different cellular processes during *L. monocytogenes* infection, at least in HeLa cells. ARPC4 was found dispensable for bacterial entry, but taking into account the central place of this subunit in Arp2/3 complex function according to previous functional and structural results [51,52], it is possible that residual ARPC4 upon RNAi treatment suffices for partial complex function.

**Vaccinia virus mobility.** In a recent study of actin polymerization by Vaccinia virus, specific roles for ARPC1B and ARPC5B have been found [10]. Indeed, it has been observed that Arp2/3 complexes containing ARPC1B and ARPC5B (named ARPC5L in this work) are significantly more efficient at promoting actin assembly than those containing ARPC1A and ARPC5A. Actin networks induced by complexes containing the subunits ARPC1B and ARPC5B were found more stable since in the presence of these specific subunits, cortactin stabilizes the Arp2/3 complexes against coronin-mediated disassembly [10].

Overall, these three reports indicate that only Arp2, Arp3 are directly partaking in the actin nucleating activity of the complex, together defining an actin nucleation core module, whereas the other subunits serve alternative roles such as determining the efficiency of actin nucleation, localization of the complex, as well as serving as an auto-inhibitory mechanism [8,51,52]. Together, the other subunits thus define the regulatory module of the Arp2/3 complex. In the case of the *L. monocytogenes* model, it is interesting to mention that vinculin inactivation by RNAi did not perturb bacterial cellular invasion nor actin-based motility [9], raising the possibility that other cellular molecule(s) not yet identified may participate to the localization/modulation of the Arp2/3 complex during *L. monocytogenes* infection-related processes.

## Concluding Remarks

While the Arp2/3 complex has been classically considered as a single molecular entity for 20 years since its discovery, an emerging possibility from recent research suggests that multiple versions of the Arp2/3 complex may co-exist in cells (**Figure 1** presents a summary of currently described complexes and their mode of regulation).



Indeed, the subunits ARPC1, ARPC3, ARPC4 and ARPC5 can be replaced by vinculin and  $\alpha$ -actinin in focal adhesions [8].

In the case of the *L. monocytogenes* system, even if dispensable, ARPC5 can be detected at both bacterial entry sites and actin comet tails by fluorescence microscopy, indicating that Arp2/3 complexes containing this subunit, while not required, may still be recruited during both processes [9]. It is possible that Arp2/3 complexes of different composition have overlapping functions during *L. monocytogenes* infection, but current data suggests that the precise composition of different Arp2/3 complexes plays a role in fine-tuning actin rearrangements in both instances. This hypothesis is supported by the observation that ARPC4 can be found predominantly early during bacterial actin comet tail formation and that knock down of ARPC4 affects initial actin polymerization at the bacterial surface rather than actin tail elongation indicating that different Arp2/3 complexes may be required in a sequential manner. The fine tuning of Arp2/3 complex actin polymerization activity depending on the subunit composition is also supported by results on the Vaccinia virus system [10].

Overall, the reports discussed herein point to the possibility that Arp2/3 is a natural modular nano-machine, capable of regulating its activity via replacement of its subunits.

## **Acknowledgments**

This work was supported by the Institut Pasteur, the Institut National de la Santé et de la Recherche Médicale (INSERM Unité 604), the Institut National de la Recherche Agronomique (INRA Unité Sous Contrat 2020), the Institut Pasteur 'Programmes Transversaux de Recherche' (PTR 460 and PTR521 to JPC), L'Agence Nationale de la Recherche (ANR-15-CE15-0017 StopBugEntry to JPC), Fondation Le Roch Les Mousquetaires, European Research Council Advanced Grant (670823 BacCellEpi to PC) and under grant agreement n° 294852-SynAd (to BG). P.C. is an International Senior Research Scholar of the Howard Hughes Medical Institute. B.G. holds the Erwin Neter Chair in Tumor and Cell Biology. The authors declare no conflict of interest.

## Figure

**Figure 1:** Diversity of Arp2/3 complexes. Central circle (gray): the canonical 7-subunit form. Top left (blue): Arp2/3 complex variants used by Vaccinia virus (alternative subunits are enclosed by a red line). Many combinations of Arp2/3 subunits are recruited by the virus. Bottom left: two “hybrid complexes”, containing the actin nucleation core and vinculin, or vinculin plus  $\alpha$ -actinin, which presumably localize the complex to focal adhesions. Bottom right: Arp2/3 complexes hijacked by *L. monocytogenes* during cellular infection (alternative subunits are enclosed by a red bold line, dispensable subunits are enclosed by a red pointed line). Top right: variations in Arp2/3 complexes caused by phosphorylation of specific subunits. The effect of the subunit substitution on the actin nucleation activity is color coded (red: reduce; green: enhance).

## Table

**Table 1:** Arp2/3 complex nomenclature

## References

- 1 Fletcher, D.A. and Mullins, R.D. (2010) Cell mechanics and the cytoskeleton. *Nature* 463, 485–492
- 2 Pollard, T.D. and Cooper, J.A. (2009) Actin, a central player in cell shape and movement. *Science* 326, 1208–1212
- 3 Berriman, M. *et al.* (2005) The genome of the African trypanosome *Trypanosoma brucei*. *Science* 309, 416–422
- 4 Goley, E.D. and Welch, M.D. (2006) The ARP2/3 complex: an actin nucleator comes of age. *Nat Rev Mol Cell Biol* 7, 713–726
- 5 Krause, M. and Gautreau, A. (2014) Steering cell migration: lamellipodium dynamics and the regulation of directional persistence. *Nat Rev Mol Cell Biol* 15, 577–590
- 6 Campellone, K.G. and Welch, M.D. (2010) A nucleator arms race: cellular control of actin assembly. *Nat Rev Mol Cell Biol* 11, 237–251
- 7 Pollard, T.D. and Beltzner, C.C. (2002) Structure and function of the Arp2/3 complex. *Curr Opin Struct Biol* 12, 768–774
- 8 Chorev, D.S. *et al.* (2014) Regulation of focal adhesion formation by a vinculin-Arp2/3 hybrid complex. *Nature Communications* 5, 1–11
- 9 Kühbacher, A. *et al.* (2015) Genome-Wide siRNA Screen Identifies Complementary Signaling Pathways Involved in Listerial Infection and Reveals Different Actin Nucleation Mechanisms during Listeria Cell Invasion and Actin Comet Tail Formation. *mBio* 6, e00598–15
- 10 Abella, J.V.G. *et al.* (2015) Isoform diversity in the Arp2/3 complex determines actin filament dynamics. *Nat Cell Biol* 18, 76–86
- 11 Machesky, L.M. *et al.* (1994) Purification of a cortical complex containing two unconventional actins from *Acanthamoeba* by affinity chromatography on profilin-agarose. *J Cell Biol* 127, 107–115
- 12 Balasubramanian, M.K. *et al.* (1996) Fission yeast Sop2p: a novel and evolutionarily conserved protein that interacts with Arp3p and modulates profilin function. *EMBO J* 15, 6426–6437
- 13 Winter, D. *et al.* (1997) The complex containing actin-related proteins Arp2 and Arp3 is required for the motility and integrity of yeast actin patches. *Curr Biol* 7, 519–529

- 14 Welch, M.D. *et al.* (1997) The human Arp2/3 complex is composed of evolutionarily conserved subunits and is localized to cellular regions of dynamic actin filament assembly. *J Cell Biol* 138, 375–384
- 15 Welch, M.D. *et al.* (1997) Actin polymerization is induced by Arp2/3 protein complex at the surface of *Listeria monocytogenes*. *Nature* 385, 265–269
- 16 Welch, M.D. *et al.* (1998) Interaction of human Arp2/3 complex and the *Listeria monocytogenes* ActA protein in actin filament nucleation. *Science* 281, 105–108
- 17 Kocks, C. *et al.* (1992) *L. monocytogenes*-induced actin assembly requires the actA gene product, a surface protein. *Cell* 68, 521–531
- 18 Yarar, D. *et al.* (1999) The Wiskott-Aldrich syndrome protein directs actin-based motility by stimulating actin nucleation with the Arp2/3 complex. *Curr Biol* 9, 555–558
- 19 Rohatgi, R. *et al.* (1999) The interaction between N-WASP and the Arp2/3 complex links Cdc42-dependent signals to actin assembly. *Cell* 97, 221–231
- 20 Machesky, L.M. *et al.* (1999) Scar, a WASP-related protein, activates nucleation of actin filaments by the Arp2/3 complex. *Proc Natl Acad Sci USA* 96, 3739–3744
- 21 Weed, S.A. *et al.* (2000) Cortactin localization to sites of actin assembly in lamellipodia requires interactions with F-actin and the Arp2/3 complex. *J Cell Biol* 151, 29–40
- 22 Mullins, R.D. *et al.* (1998) The interaction of Arp2/3 complex with actin: nucleation, high affinity pointed end capping, and formation of branching networks of filaments. *Proc Natl Acad Sci USA* 95, 6181–6186
- 23 Svitkina, T.M. and Borisy, G.G. (1999) Arp2/3 complex and actin depolymerizing factor/cofilin in dendritic organization and treadmilling of actin filament array in lamellipodia. *J Cell Biol* 145, 1009–1026
- 24 Blanchoin, L. *et al.* (2000) Direct observation of dendritic actin filament networks nucleated by Arp2/3 complex and WASP/Scar proteins. *Nature* 404, 1007–1011
- 25 Cossart, P. (2000) Actin-based motility of pathogens: the Arp2/3 complex is a central player. *Cell Microbiol* 2, 195–205
- 26 Pizarro-Cerdá, J. *et al.* (2012) Entry of *Listeria monocytogenes* in mammalian epithelial cells: an updated view. *Cold Spring Harb Perspect Med* 2,
- 27 Bierne, H. *et al.* (2001) A role for cofilin and LIM kinase in *Listeria*-induced phagocytosis. *J Cell Biol* 155, 101–112
- 28 Sousa, S. *et al.* (2007) Src, cortactin and Arp2/3 complex are required for E-cadherin-mediated internalization of *Listeria* into cells. *Cell Microbiol* 9, 2629–2643

- 29 Egile, C. *et al.* (1999) Activation of the CDC42 effector N-WASP by the Shigella flexneri IcsA protein promotes actin nucleation by Arp2/3 complex and bacterial actin-based motility. *J Cell Biol* 146, 1319–1332
- 30 Bougnères, L. *et al.* (2004) Cortactin and Crk cooperate to trigger actin polymerization during Shigella invasion of epithelial cells. *J Cell Biol* 166, 225–235
- 31 Suzuki, T. *et al.* (1998) Neural Wiskott-Aldrich syndrome protein is implicated in the actin-based motility of Shigella flexneri. *EMBO J* 17, 2767–2776
- 32 Jeng, R.L. *et al.* (2004) A Rickettsia WASP-like protein activates the Arp2/3 complex and mediates actin-based motility. *Cell Microbiol* 6, 761–769
- 33 Gouin, E. *et al.* (2004) The RickA protein of Rickettsia conorii activates the Arp2/3 complex. *Nature* 427, 457–461
- 34 Reed, S.C.O. *et al.* (2013) Rickettsia Actin-Based Motility Occurs in Distinct Phases Mediated by Different Actin Nucleators. *Curr Biol* DOI: 10.1016/j.cub.2013.11.025
- 35 Reed, S.C.O. *et al.* (2012) Rickettsia parkeri invasion of diverse host cells involves an Arp2/3 complex, WAVE complex and Rho-family GTPase-dependent pathway. *Cell Microbiol* 14, 529–545
- 36 Cudmore, S. *et al.* (1995) Actin-based motility of vaccinia virus. *Nature* 378, 636–638
- 37 Frischknecht, F. *et al.* (1999) Actin-based motility of vaccinia virus mimics receptor tyrosine kinase signalling. *Nature* 401, 926–929
- 38 Stamm, L.M. *et al.* (2003) Mycobacterium marinum Escapes from Phagosomes and Is Propelled by Actin-based Motility. *Journal of Experimental Medicine* 198, 1361–1368
- 39 Sitthidet, C. *et al.* (2010) Actin-Based Motility of Burkholderia thailandensis Requires a Central Acidic Domain of BimA That Recruits and Activates the Cellular Arp2/3 Complex. *J Bacteriol* 192, 5249–5252
- 40 Gouin, E. *et al.* (2005) Actin-based motility of intracellular pathogens. *Curr Opin Microbiol* 8, 35–45
- 41 Benanti, E.L. *et al.* (2015) Virulent Burkholderia Species Mimic Host Actin Polymerases to Drive Actin-Based Motility. *Cell* 161, 348–360
- 42 Gouin, E. *et al.* (2015) Intracellular Bacteria Find the Right Motion. *Cell* 161, 199–200

- 43 Lees-Miller, J.P. *et al.* (1992) Identification of act2, an essential gene in the fission yeast *Schizosaccharomyces pombe* that encodes a protein related to actin. *Proc Natl Acad Sci USA* 89, 80–83
- 44 Winter, D.C. *et al.* (1999) Genetic dissection of the budding yeast Arp2/3 complex: a comparison of the in vivo and structural roles of individual subunits. *Proc Natl Acad Sci USA* 96, 7288–7293
- 45 Machesky, L.M. *et al.* (1997) Mammalian actin-related protein 2/3 complex localizes to regions of lamellipodial protrusion and is composed of evolutionarily conserved proteins. *Biochem J* 328, 105–112
- 46 May, R.C. *et al.* (2000) Involvement of the Arp2/3 complex in phagocytosis mediated by FcγR or CR3. *Nat Cell Biol* 2, 246–248
- 47 Luna, A. *et al.* (2002) Regulation of protein transport from the Golgi complex to the endoplasmic reticulum by CDC42 and N-WASP. *Mol Biol Cell* 13, 866–879
- 48 Demali, K.A. *et al.* (2002) Recruitment of the Arp2/3 complex to vinculin: coupling membrane protrusion to matrix adhesion. *J Cell Biol* 159, 881–891
- 49 Kolluri, R. *et al.* (1996) Direct interaction of the Wiskott-Aldrich syndrome protein with the GTPase Cdc42. *Proc Natl Acad Sci USA* 93, 5615–5618
- 50 Kirchhausen, T. and Rosen, F.S. (1996) Disease mechanism: unravelling Wiskott-Aldrich syndrome. *Curr Biol* 6, 676–678
- 51 Gournier, H. *et al.* (2001) Reconstitution of human Arp2/3 complex reveals critical roles of individual subunits in complex structure and activity. *Mol Cell* 8, 1041–1052
- 52 Rouiller, I. *et al.* (2008) The structural basis of actin filament branching by the Arp2/3 complex. *J Cell Biol* 180, 887–895
- 53 Kelly, A.E. *et al.* (2006) Actin binding to the central domain of WASP/Scar proteins plays a critical role in the activation of the Arp2/3 complex. *J Biol Chem* 281, 10589–10597
- 54 Vadlamudi, R.K. *et al.* (2004) p41-Arc subunit of human Arp2/3 complex is a p21-activated kinase-1-interacting substrate. *EMBO Rep* 5, 154–160
- 55 LeClaire, L.L. *et al.* (2008) Phosphorylation of the Arp2/3 complex is necessary to nucleate actin filaments. *J Cell Biol* 182, 647–654
- 56 Choi, C.H. *et al.* (2013) Phosphorylation of Actin-related Protein 2 (Arp2) Is Required for Normal Development and cAMP Chemotaxis in *Dictyostelium*. *Journal of Biological Chemistry* 288, 2464–2474

- 57 Michard, C. *et al.* (2015) The Legionella Kinase LegK2 Targets the ARP2/3 Complex To Inhibit Actin Nucleation on Phagosomes and Allow Bacterial Evasion of the Late Endocytic Pathway. *mBio* 6, e00354–15
- 58 Jay, P. *et al.* (2000) ARP3beta, the gene encoding a new human actin-related protein, is alternatively spliced and predominantly expressed in brain neuronal cells. *Eur. J. Biochem.* 267, 2921–2928
- 59 Laurila, E. *et al.* (2009) Characterization of the 7q21-q22 amplicon identifies ARPC1A, a subunit of the Arp2/3 complex, as a regulator of cell migration and invasion in pancreatic cancer. *Genes Chromosom. Cancer* 48, 330–339
- 60 Millard, T.H. *et al.* (2003) Identification and characterisation of a novel human isoform of Arp2/3 complex subunit p16-ARC/ARPC5. *Cell Motil. Cytoskeleton* 54, 81–90
- 61 Boujemaa-Paterski, R. *et al.* (2001) Listeria protein ActA mimics WASp family proteins: it activates filament barbed end branching by Arp2/3 complex. *Biochemistry* 40, 11390–11404
- 62 Machesky, L.M. and Insall, R.H. (1998) Scar1 and the related Wiskott-Aldrich syndrome protein, WASP, regulate the actin cytoskeleton through the Arp2/3 complex. *Curr Biol* 8, 1347–1356
- 63 Linardopoulou, E.V. *et al.* (2007) Human Subtelomeric WASH Genes Encode a New Subclass of the WASP Family. *PLoS Genet.* 3, e237
- 64 Campellone, K.G. *et al.* (2008) WHAMM Is an Arp2/3 Complex Activator That Binds Microtubules and Functions in ER to Golgi Transport. *Cell* 134, 148–161
- 65 Zuchero, J.B. *et al.* (2009) p53-cofactor JMY is a multifunctional actin nucleation factor. *Nat Cell Biol* 11, 451–459
- 66 Rotty, J.D. *et al.* (2012) New insights into the regulation and cellular functions of the ARP2/3 complex. *Nat Rev Mol Cell Biol* 14, 7–12
- 67 Padrick, S.B. *et al.* (2011) Arp2/3 complex is bound and activated by two WASP proteins. *Proceedings of the National Academy of Sciences* 108, E472–9
- 68 Egile, C. *et al.* (2005) Mechanism of Filament Nucleation and Branch Stability Revealed by the Structure of the Arp2/3 Complex at Actin Branch Junctions. *PLoS Biol* 3, e383
- 69 Goley, E.D. *et al.* (2010) An actin-filament-binding interface on the Arp2/3 complex is critical for nucleation and branch stability. *Proceedings of the National Academy of Sciences* 107, 8159–8164



- 70 Balcer, H.I. *et al.* (2010) The p40/ARPC1 Subunit of Arp2/3 Complex Performs Multiple Essential Roles in WASp-regulated Actin Nucleation. *Journal of Biological Chemistry* 285, 8481–8491
- 71 Zaidel-Bar, R. and Geiger, B. (2010) The switchable integrin adhesome. *J Cell Sci* 123, 1385–1388